**CODMAN NEURO**

 **DePuySynthes**

**Special Bolus Needle**  
*For use with*  
**CODMAN<sup>®</sup> 3000 Series**  
**CONSTANT FLOW**  
**IMPLANTABLE PUMPS**  
**INSTRUCTIONS FOR USE**

REF AP-04013-5 and AP-04032-5

EC REP

Johnson & Johnson Medical Ltd.  
Pinewood Campus, Nine Mile Ride  
Wokingham, RG40 3EW  
United Kingdom

**TLB 3362K**

© 2002–2015 Codman & Shurtleff, Inc.

 Revised 01/15

**CE**  
**0086**

**Authorized 2003**

# ENGLISH

## IMPORTANT INFORMATION

Please Read Before Use

### Special Bolus Needle

**For use with CODMAN® 3000 Series**

### CONSTANT FLOW IMPLANTABLE PUMPS

#### INSTRUCTIONS FOR USE

REF AP-04013-5 and AP-04032-5

STERILE

EO

Rx Only

**PRECAUTION:** When performing extended infusions of more than 1 hour, through the bolus pathway of the Pump, the Pump should be emptied of the drug and filled with heparinized saline or preservative free saline.

#### Description

The Special Bolus Needle is designed for use with the CODMAN 3000 Series Pumps. These pumps are indicated for the continuous regional delivery of the following infusates with the associated route of administration indicated below:

Infusate	Indicated Route of Administration
2-Deoxy 5-Flourouridine (FUDR)	Arterial
Heparinized Saline	Arterial
Saline	Arterial/Intraspinal
Bacteriostatic Water	Arterial
Preservative Free Morphine Sulfate	Intraspinal
Baclofen	Intraspinal
Glycerin Injection	Arterial

The approved labelling for FUDR stipulates the indications, contraindications, and warnings for use of the drug in the pump.

The use of preservative-free morphine sulfate is indicated for patients with intractable chronic pain.

Bacteriostatic water or saline must be used to achieve the desired concentration of FUDR. Heparinized saline can be used during an interruption of FUDR therapy to maintain catheter patency.

Use saline to achieve the desired concentration of preservative-free morphine sulfate and during an interruption of morphine therapy to maintain catheter patency. Glycerin infusion is indicated for patients who are receiving continuous Hepatic Arterial FUDR Chemotherapy. Glycerin infusion is employed as a placebo to keep the catheter patent or to extend the refill interval for patients who require therapy interruption or withdrawal.

The CODMAN 3000 Series Pump is indicated for use in the adult population only.

#### Sterility

②

The CODMAN® Special Bolus Needle is intended for SINGLE USE ONLY; DO NOT RESTERILIZE. Use aseptic technique in all phases of handling. Codman Single Use devices have not been designed to undergo or withstand any form of alteration, such as disassembly, cleaning or resterilization, after a single patient use. Reuse can potentially compromise device performance and any usage beyond the design intent of this single-use device may result in unpredictable loss of functionality.

Codman & Shurtleff will not be responsible for any product that is resterilized, nor accept for credit or exchange any product that has been opened but not used.

As long as the inner unit is not opened or damaged, the product is sterile and nonpyrogenic.

**CAUTION:** Inspect the sterile package carefully. Do not use if:

- the package or seal appears damaged,
- contents appear damaged, or
- the expiry date has passed.

### **BOLUS PROCEDURE**

**WARNING: NEVER attempt to REFILL the Pump using a SPECIAL BOLUS NEEDLE. This use will result in giving a bolus injection to the patient and may cause a drug overdose.**



**PRECAUTION: 1) Before performing a bolus injection of any drug, review all warnings, precautions, indications and contraindications on the drug labeling.**

**PRECAUTION: 2) Use only a SPECIAL BOLUS NEEDLE for performing a bolus procedure.**

**PRECAUTION: 3) Do not aspirate fluid/blood back through the bolus path as catheter occlusion may result.**

**PRECAUTION: 4) Do not use a mechanical pressure injector system to accomplish a bolus procedure. Pressures should not exceed 40 psi when administering a bolus injection or infusion. For injections use only 10 mL (or larger) syringes and do not inject or infuse at a rate greater than 5 mL/min.**


1. Palpate the pump site and locate the raised septum.
2. Using aseptic technique, open the bolus needle pouch and expose the contents.
3. Don sterile gloves. In accordance with your institution's procedure, prep the pump site, extending the prepped area beyond the periphery of the pump. Cover the area adjacent to the pump with a fenestrated sterile drape.

4. Remove the Bolus Needle from pouch.
5. Remove cap from proximal end of set. Connect a 10 mL syringe containing saline (see CAUTION note) and flush\* set to remove air. Close clamp.

* Special Bolus Needle approximate priming volume: 0.1 mL		
Pump Bolus Pathway approximate dead volume: (Not including catheter)	CODMAN® 3000-16	0.2 mL
	CODMAN® 3000-30	0.3 mL
	CODMAN® 3000-50	0.4 mL

6. Remove needle sheath.
7. Re-palpate the pump site and locate the raised septum.

8. Insert the Special Bolus Needle into the septum until contact is made with the needle stop. Insure that the needle is PERPENDICULAR to the Pump and that contact with the needle stop is maintained throughout the procedure.

**Note:** If the needle does not remain fully inserted and in contact with the needle stop the bolus valve will close and it will not be possible to accomplish the bolus procedure.

**PRECAUTION:** When the system is flushed with saline the patient will receive a bolus dose of drug equal to the volume of drug contained in the internal bolus pathway of the Pump, plus the volume of drug in the catheter. The volume of drug in the internal pathway of the pump is shown below. The volume of drug contained in the catheter can be calculated by multiplying the length (in cm) of the catheter utilized, by the volume of the fluid contained per cm of catheter.

9. Open the clamp and flush the system with 3 mL of saline. Close the clamp and remove the syringe.

#### 10. FOR BOLUS INJECTIONS:

Attach a 10 mL syringe containing the bolus drug to the tubing set. Open the clamp and slowly inject the desired amount of drug. **Do not inject at a rate greater than 5 mL/minute.**

#### 11. FOR BOLUS INFUSIONS:

Attach a purged line from the external Pump or controller system to the SPECIAL BOLUS NEEDLE tubing set. Open the clamp on the SPECIAL BOLUS NEEDLE tubing set and begin the infusion procedure. **Do not infuse at a rate greater than 5 mL/minute.**

**PRECAUTION:** For infusions of more than 1 hour, use the SPECIAL BOLUS WINGED INFUSION SET AP-04033 or AP-04034.

**PRECAUTION:** The patient should remain still throughout this procedure to prevent any movement of the SPECIAL BOLUS NEEDLE. If the needle is withdrawn during the procedure, the bolus safety valve will close and it will not be possible to accomplish this procedure. If the needle is withdrawn an amount exposing the lumen above the upper septum of the Pump, a drug extravasation may occur.

12. After the drug injection or infusion, close the clamp, remove the drug syringe or line from the external Pump and replace it with the 10 mL syringe containing saline. Open the clamp and slowly inject 3 mL of saline.
13. Close the clamp, remove the needle from the patient and apply pressure and a bandage to the puncture site.

- CAUTION:**
- A) Always use preservative free saline for flush solution with pumps that access the intrathecal or intraspinal spaces.
  - B) If more than one drug is administered, flush set and pump with saline between drugs.
  - C) SPECIAL BOLUS NEEDLE is not intended for use with Implantable Ports.

- D) SPECIAL BOLUS NEEDLE is intended for single use only. DO NOT RESTERILIZE.
- E) As with all conventional (PVC) administration sets, this device should be used in compliance with the drug/solution labeling.

.....

**Straight Special Bolus Needles**

Specification	REF No.	Packaged
20 gauge, 1.4 inch (36 mm)	AP-04013-5	5 units/case
20 gauge, 2.0 inch (51 mm)	AP-04032-5	5 units/case

© CODMAN is a registered trademark of Codman & Shurtleff, Inc.

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