Codman

CODMAN®
Pump Catheter
(REF 60-2914) for
MedStream™ and
ARCHIMEDES®
Pumps

EC REP

Codman

A division of Johnson & Johnson Medical Ltd.
Pinewood Campus, Nine Mile Ride
Wokingham RG40 3EW
United Kingdom

CE
0086

Authorized 2003
LCN 198830-001/J

© 2002–2010 Codman & Shurtleff, Inc.

Recycled Revised 07/10

ENGLISH

IMPORTANT INFORMATION

Please Read Before Use

CODMAN® Pump Catheter (REF 60-2914) for MedStream™ and ARCHIMEDES® Pumps

STERILE	EO
---------	----

Refer to the user's manual provided with the infusion pumps for additional information.

Indications

The CODMAN Pump Catheter is indicated for use as part of an infusion system using a MedStream or ARCHIMEDES Infusion Pump. The pumps are indicated for the continuous intrathecal infusion of selected drugs for pain management or relief of spasticity.

Contraindications

Active or incompletely treated infection.

WARNINGS

The use of non-Codman accessories with the pump catheter can result in damage to the catheter or failure of the system to function as intended.

Improper use of implanted, constant-flow infusion pumps could result in drug under- or overdose. Users must comply with the product instructions for initial pump preparation, implantation, initial filling, refilling, and injecting into the bolus port of the pump.

Only qualified medical personnel must perform these tasks.

Before closing the incisions, verify that the central port and bolus port of the implanted pump can be palpated, that the catheter(s) will not become twisted or contorted, and that the catheter(s) will not obstruct the central port or the bolus port.

Never inject into the pump catheter; a drug under- or overdose or damage to the catheter can result.

Precautions

Inspect the sterile package carefully. Do not use if:

- the package or seal appears damaged;
- the contents appear damaged; or
- the expiry date has passed.

The pump catheter is **for single use only. Do not reuse.**

Use sterile technique in all phases of handling this product.

Exercise extreme care to prevent the silicone catheter(s) from coming in contact with towels, drapes, talc, or any linty or granular surfaces. Silicone is highly electrostatic and as a result, attracts airborne particles and surface contaminants that could produce tissue reactions.

Silicone has a low cut and tear resistance; therefore, do not tie ligatures too tightly on the strain relief sleeves or the connectors. Do not use stainless steel or chromium ligatures. Use only No. 0 silk suture for securing the catheters.

Do not handle catheters with sharp instruments. Use rubber shod forceps when handling catheters. Take care not to inadvertently cut or puncture catheters.

Make sure that catheter placement and connections are secure. Failure to adequately connect and secure catheters in place can result in:

- dislodgment or disconnection of the catheter;
- obstruction of the central port or bolus port;
- cessation of therapy; or
- delivery of drug to the pump pocket or the subcutaneous tissue.

It is important to note the total catheter length (pump catheter and/or intraspinal catheter) and the catheter inner diameter in the patient's record. The total implanted catheter length and

inner diameter are used to calculate the time required for drug to advance to the intraspinal catheter tip, as well as to help prevent a drug overdose when injecting into the bolus port. Refer to *Catheter Volume Calculation*.

Do not trim catheter(s) before placement.

If a catheter is to be trimmed, leave an extra length of catheter to help ensure that no tension will be placed on the catheter. Measure and record the length of catheter removed. This information is necessary for calculating the volume of the implanted catheter.

During implantation, verify that catheter(s) will not become kinked or occluded due to knots, tight geometries, or a tortuous position. To prevent dislodgment or kinking, always use an anchor to secure the intraspinal catheter to the surrounding tissue.

Pass each catheter before connecting whenever possible. If passing a catheter that is already connected or implanted, ensure that the direction of tunneling and passing the catheter leads away from the site of connection or implantation.

Use only No. 0 silk suture for securing catheters and connectors. Smaller diameter suture and stainless steel suture may cut or tear silicone catheters.

Adverse Events

Catheter performance problems, such as kinking, disconnection, leakage, breakage, occlusion (complete or partial), dislodgment, migration, formation of fibrosis or hygroma, which can cause:

- delivery of the drug to an inappropriate site, such as the pump pocket or subcutaneous tissue
- drug withdrawal symptoms
- underinfusion of the drug
- a return of underlying symptoms
- free-floating catheter in the cerebrospinal fluid (CSF)
- CSF leakage, subcutaneous collection of CSF, or central nervous system (CNS) pressure changes
- damage to the spinal cord

- hemorrhage
- the need for surgical replacement of the catheter
- organ failure
- stroke
- death

Product Description

The CODMAN Pump Catheter is designed only for use in a system that includes a MedStream or an ARCHIMEDES Infusion Pump. The pump catheter provides the pathway from the pump to the intraspinal catheter and/or the trial catheter. The pump connector end of the pump catheter fits onto the outlet of the pump.

The pump catheter is made of medical-grade silicone elastomer impregnated with barium sulfate. **Note:** Although the catheter is X-ray detectable, on X ray a gap of 1–2 mm appears between the catheter and the pump body. This gap does not represent a faulty connection. The catheter dimensions are as follows:

Inner diameter:	0.6 mm
Outer diameter:	2.4 mm
Length:	100 cm

The pump catheter is compatible with the following catheters and connectors:

Table 1 Compatible Catheters and Connectors	
Catheter	Connector
Silicone Catheter System	
CODMAN® Intraspinal Catheter (REF 60-2918)	CODMAN Intraspinal Catheter Connector with Strain Relief Sleeve (REF 60-2950)
CODMAN Trial Catheter (REF 60-2923)	
Titanium-reinforced Catheter System	
SureStream™ Intraspinal Catheter (REF 70020)	SureStream Replacement Connector with Strain Relief Sleeve (REF 70029)

Sterility

2

This product is intended for **SINGLE USE ONLY; DO NOT RESTERILIZE.**

Codman Neuro Sciences Sàrl will not be responsible for any product that is resterilized, nor accept for credit or exchange any product that has been opened but not used.

As long as the inner unit is not opened or damaged, the product is sterile and nonpyrogenic.

Instructions for Use


A. Attaching the pump catheter to the pump before implantation


Refer to the pump instruction manual to perform the pump preparation. When using a MedStream pump, also refer to the MedStream Infusion System Programming Guide.

Refer to the instructions provided with the intraspinal catheter to implant it and join it with the pump catheter.

1. Remove the protective tubing from the outlet of the pump, if applicable.
2. Gently push the pump connector end of the catheter onto the pump outlet with a twisting motion.
3. Make sure the groove of the connector is completely seated in the groove of the pump outlet.
4. Use No. 0 silk suture material to tie the catheter securely to the pump outlet. The suture must lie in the groove of the connector. See Figure 1.

Fig. 1

- A. Pump outlet
- B. Needle stop
- C. Catheter connector
- D. MedStream pump
- E. Catheter
- F. Ligature

5. Follow the user manual packaged with the pump to complete the pump preparation, functional test, and implantation.
6. Use the catheter connector and, if appropriate, the strain relief sleeve to join the pump catheter to the intraspinal catheter or to the trial catheter. Follow the instructions for use packaged with the catheter connector and strain relief sleeve or with the intraspinal catheter kit.

B. Replacing an implanted pump catheter

Required materials:

- Appropriate Bolus Kit:
 - For MedStream pump** REF 91-4284
 - For ARCHIMEDES pump** REF 60-2922
- No. 0 silk suture
- 10 mL preservative-free sterile 0.9% saline solution
- Appropriate replacement catheter connector (refer to Table 1)

1. **For MedStream pump only.** Use the MedStream Control Unit to interrogate the MedStream pump. (Refer to the MedStream Infusion System Programming Guide for more information.) Scroll to the catheter information section (the second page of the summary screens). Make note of the “Total Catheter Volume” of the current catheters.
2. Referring to the instructions provided with the bolus kit, use the bolus needle to aspirate the contents of the catheters.

3. Clamp the intraspinal or trial catheter near its connection with the implanted pump catheter.
4. Cut the intraspinal or trial catheter just distal to the catheter connector. Record the length of catheter removed.
5. Remove and discard the old pump catheter.
6. Gently push the pump connector end of the new pump catheter onto the pump outlet with a twisting motion.
7. Make sure the groove of the connector is completely seated in the groove of the pump outlet.
8. Use No. 0 silk suture material to tie the catheter securely to the pump outlet. The suture must lie in the groove of the connector. See Figure 1.
9. To test patency of the catheter and connection and to prime the pump catheter, administer a 5 mL bolus injection of saline through the bolus port. Refer to the instructions for administering a bolus in the instructions provided with the bolus kit. Resistance may be a sign of occlusion; take appropriate measures. Inspect the pump connection for signs of leakage.
10. Trim the pump catheter as appropriate. Use the appropriate catheter connector and strain relief sleeve to join the pump catheter to the intraspinal catheter. Follow the instructions packaged with the catheter connector. Record the amount of pump catheter removed.
11. Remove the occluding clamp from the intraspinal catheter.
12. Administer another 5 mL bolus injection of saline to test the connection of the catheters. Resistance may be a sign of occlusion; take appropriate measures. Inspect the connection of the catheters for signs of leakage.

- 
13. **For MedStream pump only:** use the Control Unit to update the catheter length information. Refer to the Programming Guide for additional information.
 14. If desired, deliver a bolus of drug solution equal to the volume of both catheters to clear the saline and begin immediate intraspinal delivery of drug therapy.

WARNING: This revision procedure results in the catheters being filled with saline solution. Intraspinal delivery of the drug solution from the drug reservoir will be delayed until the saline solution exits the intraspinal catheter. Carefully calculate the amount of the delay, based on the total catheter length and the flow rate, before administering additional medication by another delivery method.

For MedStream pump only: use the steps below and refer to the MedStream System Programming Guide to program a bolus.

- 
- 
- a. Interrogate the pump and click on “Catheter Info” at the bottom of the pump information screen to read the catheter information. Make note of the volume of both catheters.
 - b. Calculate the bolus dose by using the formula:
Bolus dose = (pump catheter volume + intraspinal catheter volume) x drug concentration
 - c. Follow the instructions in the “Single Bolus” section of the Programming Guide to program a bolus.

Catheter Volume Calculation

Length of catheter (cm)	Volume (mL) (0.6 mm Inner Diameter)
20	0.057
30	0.085
40	0.113
50	0.141
60	0.170
70	0.198
80	0.226
90	0.254
100	0.283

Warranty

Codman & Shurtleff, Inc. warrants that this medical device is free from defects in both materials and workmanship. **Any other express or implied warranties, including warranties of merchantability or fitness, are hereby disclaimed. Suitability for use of this medical device for any particular surgical procedure should be determined by the user in conformance with the manufacturer's instructions for use. There are no warranties that extend beyond the description on the face hereof.**

® CODMAN and ARCHIMEDES are registered trademarks of Codman & Shurtleff, Inc.
™ MedStream and SureStream are trademarks of Codman & Shurtleff, Inc.

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