# Codman

## Model 400 and 550 Implantable Pump

### Refill Kit

REF IP-40820

EC REP

## Codman

A division of Johnson & Johnson Medical Ltd.  
Pinewood Campus, Nine Mile Ride  
Wokingham RG40 3EW  
United Kingdom

CE  
0086

Authorized 2003

TLB 3250/E

© 2002–2013 Codman & Shurtleff, Inc.

Revised 02/13

# ENGLISH

## IMPORTANT INFORMATION

Please Read Before Use

### Model 400 and 550 Implantable Pump Refill Kit

REF IP-40820

STERILE EO

Rx Only

#### The Refill Kit contains:

- 1 – calibrated syringe barrel
- 1 – tubing set with stopcock
- 4 – (22 gauge)\* Huber Point Needles
- 3 – alcohol swabsticks
- 3 – iodine swabsticks
- 6 – gauze pads
- 1 – fenestrated drape
- 2 – pump templates
- 2 – adhesive bandages

\*Two different lengths (1" and 1 1/2") Huber Point Needles are provided. For pump septum access, choose the most suitable length.

Not provided in kit:

- 1 – pair of sterile gloves
- 1 – syringe of appropriate infusate, warmed to room temperature (use smallest available syringe size which will accommodate refill volume):  
MODEL 400: 50 ml (single and dual)  
MODEL 550: 32 ml

#### How Supplied

②

The CODMAN® Model 400 and 550 Implantable Pump Refill Kit is intended for SINGLE USE ONLY; DO NOT RESTERILIZE. Use aseptic technique in all phases of handling. Codman Single Use devices have not been designed to undergo or withstand any form of alteration, such as disassembly, cleaning or resterilization, after a single patient use. Reuse can potentially compromise device performance and any usage beyond the design intent of this single use device may result in unpredictable loss of functionality.

Codman & Shurtleff will not be responsible for any product that is resterilized, nor accept for credit or exchange any product that has been opened but not used.

As long as the individual package is not opened or damaged, the product is sterile.

Testing has shown that the following components are nonpyrogenic.

- Tubing set
- Stopcock
- Caps
- Syringe
- Needles

**CAUTION:** Inspect the sterile package carefully. Do not use if:

- the package or seal appears damaged,
- contents appear damaged, or
- the expiry date has passed.

#### A. Refill Procedure


##### CAUTIONS

- The health care team should be thoroughly familiar with the pump refill procedure and precautions before proceeding.
- It is critical to the safety and function of the pump that only 22 gauge Huber Point Needles be used. They are intended for SINGLE USE ONLY. Needle sizes larger than 22 gauge may compromise septum integrity.


- DO NOT ASPIRATE FLUID FROM THE PUMP. ASPIRATION MAY CAUSE PUMP FAILURE.
- Verify that the amount of infusate in the refill syringe is compatible with the pump reservoir volume.
- Throughout this procedure, the Clinician should confirm that all of the connections are secure in order to eliminate any chance of leakage.
- The Clinician should be certain that the needle is properly positioned in the center septum of the pump prior to any drug instillation.
- It is essential that accurate refill records be maintained.

**NOTE:** Refer to the Clinician's Manual for appropriate drug and flushing solutions.


1. Using sterile technique, expose the components of the Refill Kit to the sterile field.
2. Identify the outer perimeter of the pump by palpating the pump pocket.
3. Don sterile gloves. Use alcohol swabsticks to prep the pump site in a circular fashion, beginning at the center septum and working outward. Extend the prepped area beyond the periphery of the pump. Use iodine swabsticks and repeat.

4. Place the fenestrated drape over the prepped pump site.
5. Position the appropriate pump template over the pump (optional).
6. Attach the empty syringe barrel to the stopcock and a 22 gauge Huber Point Needle. TIGHTEN ALL CONNECTIONS. Place on sterile field.


7. Attach a second 22 gauge Huber Point Needle to the luer lock connection of the tubing. Insert this needle into the center septum of the pump at a perpendicular angle.

**NOTE:** Use of the tubing injection site is optional. To remove, unscrew the luer lock which adjoins the tubing. Stopcock and

empty syringe barrel should then be attached directly to tubing set prior to penetrating pump septum (Skip step #8).

8. Insert the needle with stopcock and empty syringe barrel attached into the injection site of the tubing.
9. With stopcock in the open position, allow the pump to completely empty. Close the stopcock. Note the amount of returned volume in the syringe barrel (the syringe barrel has been precalibrated to include tubing volume).


**NOTE:** If no fluid returns into the syringe barrel when the center septum is accessed, the possibility exists that:

- a. the needles or tubing set may be blocked;
  - b. the pump is empty of infusate; or
  - c. the septum may not have been penetrated.
- Leaving the needle in the center septum, disconnect the tubing set. Remove the syringe barrel and stopcock and place on the sterile field. Attach a 10ml syringe of appropriate placebo to the needle in the tubing injection site. (If tubing injection sites not used, leave stopcock on tubing and attach placebo syringe to stopcock.) Prime tubing, leaving 5ml of solution in the syringe. Reconnect the tubing set to the needle in the center septum and inject the 5ml into the pump. Release the pressure on the syringe plunger. If fluid returns into the syringe, indications are that the pump was empty prior to refill.

If no fluid returns into the syringe, reattempt to locate and penetrate the septum. Repeat the 5ml injection of placebo.


If repeated attempts to locate the septum do not result in a fluid return into the syringe, contact Codman Customer Service Department for assistance.

10. Close clamp. Remove the syringe barrel, stopcock, and needle from the injection site of the tubing and set aside. **NOTE:** If tubing injection site is not used, remove syringe barrel and stopcock from tubing.


11. Attach a needle to the syringe of infusate and insert into the tubing injection site. **NOTE:** If tubing injection site is not used, attach syringe of infusate directly to tubing set.

- Open clamp. Using both hands, inject 5ml of infusate into the pump. Release pressure on the plunger and allow a small amount of fluid to return to the syringe. This test confirms proper needle placement. Continue to inject and check needle placement at 5ml increments until the syringe is emptied.

- Maintaining positive pressure on the plunger, clamp the tubing as close as possible to the needle in the center septum.

- Gently pull the needle out of the center septum and apply pressure to the puncture site. Remove iodine — apply adhesive bandage.

## B. SIDEPORT™ Bolus Procedure

- Constant Flow Pumps

The Refill Kit contains supplies required for this procedure.

Additional supplies needed:

- 2 – 10ml syringes containing appropriate flushing solution
- 1 – drug syringe (limit syringe size to 10ml or larger)
- 1 – FOR INTRA-ARTERIAL APPLICATIONS ONLY: 10ml syringe filled with 3ml of heparinized normal saline (100U/ml)

### CAUTIONS


- Only 22 gauge Huber Point Needles can be used with Models 400 and 550 Implantable Pumps. These needles are intended for SINGLE USE ONLY. Needle sizes larger than 22 gauge may compromise septum integrity.
- The health care team should be thoroughly familiar with the pump SIDEPORT bolus procedure and precautions before proceeding.
- Prior to injection into the SIDEPORT, review all warnings, precautionary notes, indications and contraindications for radiolabeled and therapeutic drugs.
- Misdirected needle placement during SIDEPORT access may damage the segment of the catheter in the pump pocket.

- POSITIVE PRESSURE MUST ALWAYS BE MAINTAINED ON THE SYRINGE PLUNGER WHEN ACCESSING THE SIDEPORT.
- DO NOT ASPIRATE FROM THE SIDEPORT(S) AT ANY TIME. ASPIRATION MAY CAUSE PERMANENT OCCLUSION OF THE DRUG FLOW PATH NECESSITATING EXPLANATION OF THE DEVICE.
- The injection rate must be limited to 10ml/minute and the syringe size to 10ml or larger to avoid an overpressured condition occurring in the SIDEPORT. Overpressurization could cause permanent catheter damage.
- Bolus administration through the SIDEPORT will result in the patient receiving the volume of infusate residing in the catheter.
- Confirm compatibility of bolus solution with the drug in the center reservoir of the pump.

**NOTE:** Refer to the Clinician's Manual for appropriate flushing solutions.

- Fill two 10ml syringes with an appropriate solution to flush catheter. Fill drug syringe (limit syringe size to 10ml or larger).
- Using sterile technique, expose components of the Refill Kit to the sterile field.
- Identify SIDEPORT septum by palpating the outer perimeter of the SIDEPORT.
- Don sterile gloves. Use alcohol swabsticks to prep the SIDEPORT site in a circular fashion extending the prepped area beyond the periphery of the SIDEPORT. Use iodine swabsticks and repeat.
- Place the fenestrated drape over the prepped site.
- Position the appropriate pump template over the pump (optional).
- Attach a 10ml syringe containing appropriate flushing solution to a 22 gauge Huber Point Needle. Penetrate the injection site of the tubing set. Flush through the tubing to remove any air. CLAMP THE TUBING.

**NOTE:** Use of the tubing injection site is optional. To remove, unscrew the luer lock which adjoins the tubing. The syringe and stopcock should then be attached directly to the tubing set. After the tubing is flushed, close the stopcock and clamp tubing.


- Insert the 22 gauge Huber Point Needle on the end of the tubing into the SIDEPORT septum at a perpendicular angle. MAINTAINING POSITIVE PRESSURE ON THE SYRINGE, OPEN THE TUBING CLAMP (AND STOPCOCK). Flush the pump catheter with the contents of the 10ml syringe containing placebo.

### CAUTIONS

- Incorrect needle placement may cause catheter puncture.
- Failure to maintain positive pressure could result in permanent pump or catheter occlusion necessitating explantation of the pump.

**IMPORTANT:** Confirm needle placement by examining the puncture site closely for any subcutaneous infiltration.

- While maintaining pressure on the syringe plunger, CLOSE THE CLAMP. Remove the syringe and needle from the tubing injection site.

**NOTE:** If tubing injection site is not used, close clamp and stopcock; remove syringe.

- Attach the drug syringe to a 22 gauge Huber Point Needle. Penetrate tubing injection site.

**NOTE:** If tubing injection site is not used, attach drug syringe to stopcock. Open stopcock.

While maintaining pressure on the syringe plunger, OPEN THE CLAMP. Slowly inject the contents of the drug syringe.


**CAUTION:** Injection rate must be limited to 10ml/min. to avoid an overpressurized condition. Overpressurization could cause permanent catheter damage.

- While maintaining pressure on the syringe plunger, CLOSE THE CLAMP. Carefully remove the drug syringe and needle from the tubing injection site to avoid drug spillage.

**NOTE:** If tubing injection site is not used, close clamp and stopcock; remove syringe.

- Attach the second 10ml syringe containing flushing solution. Maintaining pressure on the syringe plunger, OPEN THE CLAMP (and stopcock) and flush the catheter. Continuing to apply pressure to the syringe plunger, CLOSE THE TUBING CLAMP. Gently withdraw the needle from the SIDEPORT septum. Remove iodine—apply adhesive bandage.

**IMPORTANT:** For intra-arterial applications, catheter should be flushed with 3ml of heparinized normal saline (100U/ml) prior to withdrawing needle from septum.

® CODMAN is a registered trademark of Codman & Shurtleff, Inc.

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